



New Product!

Outstanding Support, Stability and Comfort. That's the StrongArm difference...

"From our first glimpse of the StrongArm, we knew that the concept and design fulfilled an unmet consumer need. We've had a tremendous positive consumer and therapist response when we've demonstrated the StrongArm. Comments from people using the StrongArm have been extremely favorable and requests for how to purchase the product have been continuous."

James Leahy, University at Buffalo SUNY/Department of Rehabilitation Science



StrongArm Mobility Order Numbers

StrongArm Forearm Cane - Standard
 Fits Adults up to 6'1"
 3.125" Cradle Opening
 300 lbs. User Capacity

#40011 Standard Silver Matte, 6/CS 

#40014 Standard Copper, 6/CS 

StrongArm Forearm Cane - Stout
 Fits Adults up to 6'4"
 3.5" Cradle Opening
 300 lbs. User Capacity

#50017 Stout Silver, 6/CS 

To find your closest dealer or to order a StrongArm directly, please contact us by:

Phone: 888 482 1569

Fax: 847 735 9290

or on the web at:

www.StrongArmMobility.com

The StrongArm name and logo are registered trademarks of StrongArm Mobility. Products are protected by US and foreign patents. Manufacturer reserves the right to change specifications without notice and obligation. StrongArm Forearm Canes are designed and manufactured by StrongArm Mobility, Evanston, IL USA (www.StrongArmMobility.com)



DESIGN

Features

- Cushioned forearm cradle stabilizes wrist
- One-piece upper tube
- Centered cradle
- Small radius bends
- Off-set handle design
- Reinforced tubular construction
- Aluminum parts made from high strength alloy

Benefits

- StrongArm feels like solid extension of the arm
- No rattling assemblies or confining plastic cuffs to deal with
- Insures neutral wrist position for comfort in either hand
- Elegant more compact design
- Centers weight over cane reducing stress on the elbow
- Light weight doesn't fatigue user. Solid feel, even at upper weight limit
- Ensures stability and long product life

100% Latex Free

Support

Dual Cushioned Surfaces


Sueded, latex-free foam cushions provides all-day comfort and assist user rising from chair or bed



Offset Bend

Centers grip over tip for maximum stability when cane is at user's side

Comfort



Hexagonal Floor Tip

Angled sides "flatten" to form larger footprint for improved traction on all surfaces

Forearm Bracing Cradle

Tubular "cradle" braces forearm on 3 sides to fully stabilize wrist



Ergonomic Hand Grip

Contoured grip locates hand for reduced fatigue

Stability

Improved Tightening Nut

Large spaced "lugs" around perimeter are easier to grasp



High Strength Aluminum Components

Provide lightweight strength and attractive anodized finishes

Medicare reimbursable under HCPCS codes E0111 (single), or E0110 (in pairs)

To learn more about the StrongArm Forearm Cane visit us at

StrongArmMobility.com